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

Technical information

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압전소자 특성 Characteristics of piezo-electric element

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
A piezo-electric element is a solid actuator subject to displacement when voltage is applied. Applying voltage with electrodes formed on both faces of the piezo-electric ceramic plate will cause distortion by about 0.1%.

Figure 1 shows the structure of stacked piezo-electric elements. Form many internal electrodes and stack them on the surface opposed to the piezo-electric ceramic thin plate.

On the sides, alternately fit an insulator to every two layers.

The piezo-electric element sides are covered with resin for protection. Some are packed in metal sealed case that shuts off outside air.

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Features of stacked piezo-electric element include:

- ① large generative force; ② high resolution; ③ high response; and ④ high efficiency of electric-mechanical conversion.

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Figure 2 shows the displacement caused when voltage is applied to stacked piezo-electric element. Although the piezo-electric element expands and contracts according to the voltage, the displacement curve with increasing voltage (pressurizing) and that with decreasing voltage (depressurizing) are not on the same track. In this phenomenon, which is called hysteresis, about displacement of 15 to 20% of the maximum displacement occurs.

In addition, when looking at the displacement curve after repeated application of voltage onto the piezo-electric element, you will see that the displacement in proximity of voltage 0 V varies between the first one and subsequent ones. This phenomenon is called the zero shift.

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Figure 3 shows the temporal change of displacement when retained with certain voltage applied onto the piezo-electric element. You may see that the displacement is gradually changing over time. This phenomenon is called creeping, which is caused by polarization of piezo-electric ceramics.

The creeping amount becomes larger according to the difference of input voltage from the previous voltage value and the direction will be the same as increase/decrease of the voltage.

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Due to the characteristics such as hysteresis and creeping, displacement of the piezo-electric element cannot be determined only in accordance with the applied voltage value. When the piezo-electric element is used as a driving source of a precision positioning device, it is necessary to also use a displacement sensor to measure the actual displacement (travel distance) while always controlling the applied voltage.

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A force that applies voltage to restore the original position is called the maximum generative force.

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

Displacement of the piezo-electric element varies depending on whether the applied pressure (load) is a constant load like a weight or a variable load that changes gradually like a spring.

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


If the load (pressure) is constant, the rated displacement can be achieved by the rated voltage when the contracted position is taken as a reference point.

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If the load (pressure) is an elastic force like a spring force, displacement only occurs in portions where rigidity of the piezo-electric element and that of the spring are balanced, which results in a displacement disorder.

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The type enclosed in a metal case is designed to be less susceptible to the atmosphere by having a structure in which the piezoelectric element is isolated from the external air. It is more durable to the ambient temperature and humidity than the piezoelectric element incorporated in the type with a resin exterior. In addition, this type is less subject to deterioration (life shortening) even if it is applied with a constant voltage for a long time, and is therefore suitable for the following equipment, environments or usages.

- Equipment or microscope that handles liquids such as bio
- Environment with oil or mist due to processing machinery
- Cases where the stage is used in static (stand-by) state for a long time in a certain position
* The operating time is shorter.
- Use in a place where it is difficult to replace (remove) the stage in the case of failure
- Applications where high reliability is required
- High / ultrahigh vacuum

압전소자의 절연저하 Decreased insulation of the piezoelectric element

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
One of the failure modes of a laminated piezoelectric element is short-circuiting mode (dielectric breakdown) due to reduced insulation. The mechanism of such an insulation reduction has not been completely elucidated. However, it has been confirmed that the failure incidence ratio is substantially different between static use (direct current is applied) and dynamic use (pulsing voltage is applied). It has also been confirmed that a laminated piezoelectric element is greatly affected by the humidity in addition to the applied voltage and the environmental temperature, as with general electronic parts. Moisture absorption in the element due to the ambient environment during storage or operation may cause the reduction of the insulation resistance resulting from migration*.

* Migration

It is a phenomenon in which Ag in the internal electrode is ionized, and the ionized Ag penetrates into the element along the grain boundary due to the applied electric field, resulting in dielectric breakdown. It tends to occur in cases where a constant voltage is continuously applied in a high-humidity environment.

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
For nanometer order positioning, effects of backlash and friction that are normally negligible becomes relatively significant and not negligible. Therefore, a guidance mechanism and magnification mechanism that utilize its elastic deformation are use in Piezo Stage. For example, a piezo stage with the elastic hinge as indicated in the above figure may maintain piezo's high resolution and also realizes a large stroke, ultra-small tilt and high straightness at the same time by using its magnification mechanism and spring-loaded guidance mechanism.

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Example of an actual elastic hinge mechanism and stress analysis

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- With no frictional wear, the running accuracy will not decrease even in long-term repetitive reciprocating motion.
- Since no grease is used, the product will not generate dust, and can be used on a maintenance-free basis.
- The hinge design will be conducted while considering the repeated fatigue strength through stress analysis based on FEM analysis.

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
In the catalog standard product, aluminum alloy or steel is used as the material of the body. However, a product fabricated with material that suits the application can also be available.

- Aluminum alloy : Lightweight and well workable, and frequently used also at our company.
- Steel : Being highly rigid, suitable for applications where a load bearing capacity is required.
- Stainless steel : Mainly martensitic stainless steel is used. Suitable for environment-resistant and non-magnetic applications.
- Titanium alloy : Used if the product is to be incorporated in equipment that requires non-magnetic property.
- Super invar : Since its thermal expansion coefficient is low, this material is used in cases where the impact of change in the ambient temperature needs to be minimal.

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FEM analysis of the thermal strain of aluminum alloy and super invar (an example where stages with an opening are compared)

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A capacitance type displacement sensor is a displacement sensor that measures the changes in capacitance between the probe and target. This sensor is featured by high resolution, high reproducibility and high stability.

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
The sensor output becomes zero when the distance between the sensor probe and the measurement target is a certain distance. The distance is called "offset." Set the measurement range based on the offset distance. Offset is normally equivalent to 5 to 10% of the measurement range and the measurement range that can be selected varies depending on external diameter of the measurement probe.

Sensor's displacement detection limit (resolution) is calculated by the following formula.

$$\text{Resolution} = \text{Measurement Range} / \text{Output Voltage Width} \times \text{Noise Level}$$

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The capacitance type displacement sensor adjusts the gain and linearity separately by the combination of the probe and amplifier. Therefore, it is necessary to use a stage with an integral probe and Piezo Controller with an amplifier in a one-to-one pair. Wrong combination may cause a failure or impair the performance.

Especially when a multiple-axis stage is used, channels connected to the controller are predetermined for each axis, so care must be exercised to connect in right combination.

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The open loop operations simply amplify the command signal 0-10 V fifteenfold and apply voltage of 0-150 V to the Piezo Stage onto the Piezo Stage for driving. Therefore, characteristics of piezo-electric element such as hysteresis and drift take effect as properties of the stage, resulting in positioning errors. In closed loop operations, the Piezo Controller always monitors signals from the built-in displacement in the Piezo Stage and compensates and corrects them so that the command signals 0-10 V become consistent with the stage stroke 0 to maximum μm .

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

[Example] Described below is how open loop and closed loop operations act when analog command voltage is linearly applied to the piezo stage of the following specifications in order of "0 V > 10 V > 0 V."

Travel range : 120 μm (open loop)
 100 μm (closed loop)
 Command voltage : 0→10 V→0 V

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
The standard length of "piezo cable" and "sensor cable" is 2 m. However, for cable run within the device or increase of distance between the piezo Stage and Controller, it may be necessary to extend the cable length. At this point, extending the cable will increase the electric noise, which affects the resolution and repetitive positioning accuracy. Especially the cable length affects electric signals of the capacitance displacement sensor as they are is subtle by their nature.

Therefore, the sensor amplifier board can be apart from the controller and serve as a relay sensor amplifier to amplify the sensor signals near the piezo Stage. After passing the relay sensor amplifier, the signals may be connected to the controller via a cable of any length. This allows for extension of the cable length while maintaining the original performance of the piezo Stage without impairing the S/N ratio.

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

Although products indicated in "General Catalog" are assumed to be used in the atmosphere, they can also be used in vacuum with some vacuum specification parts.

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
For vacuum specifications, mainly the following changes must be implemented.

- ① A vacuum-type piezo cable line shall be used.
- ② A vacuum-type sensor cable line shall be used.
- ③ The vacuum cable length and atmosphere cable length shall be determined.
(Ex) 0.5 m in vacuum and 1.5 m in atmosphere
- ④ BNC connectors shall be used for the ends of the piezo cable and sensor cable.
- ⑤ Adhesive for vacuum devices shall be used.
- ⑥ The surface treatment method for the Piezo Stage main unit shall be determined

Work(s) by customer

- ⑦ Please use "insulated BNC (female) terminals" as vacuum current feedthrough for the vacuum partition.
(2 BNC terminals are required for 1 axis).

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Stacked piezo-electric elements are mechanically characterized by strong compression force and weak tensile and shear strength due to alternately-stacked layers of thin ceramic and electrodes. Although they are designed so that certain preload is applied to them when being assembled onto the Piezo Stage, care must be exercised not to give any external force in pull or shear directions in actual handling of the Piezo Stage.

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

In dimensional drawings of "Piezo Stage" provided in this catalog, a direction of the trolley table surface travel is indicated by an arrow. (See the following figure.) The arrow not only represents the travel direction of the table surface under voltage but also indicates the direction of piezo-electric element expansion. That is, applying an external force to the direction of the arrow also gives a pull force on the piezo-electric element, which may result in any malfunction.

Especially, care must be exercised so that the "torque" and "torsion" generated when a work is fixed to the trolley table by screws do not pull the piezo-electric elements.

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
When installation or release, power shall not be given in the moving direction

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Care must be exercised not to make the "direction of gravitational force" a pull force for the piezo-electric elements even when the unit is used in other than the horizontal posture.

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
Please install the unit so that the Piezo Stage moves upward (under voltage). The same shall apply to installation on an inclined plane. There are no constraints in installation on a wall if it travels horizontally in parallel with the floor surface.

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The Z-axis stage cannot be hanged upside down. Customization is required for such usage.

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How to read the specification table of piezo stage

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
- ① Travel range
Maximum travel distance for closed loop and open loop.
*For the open loop, the value is a representative value.
- ② Resolution
Resolution for the closed loop control.
- ③ Repeatability
Repetitive positioning accuracy for the closed loop control.
- ④ Load capacity
The normal operating carry load is indicated as a representative value.
- ⑤ Resonant frequency
The mechanical resonance frequency under a prescribed carry load is indicated as a representative value.
- ⑥ Pitching/Yawing/Rolling
Maximum angular difference in the entire travelling range possible with the open loop.
(For 2 or 3 axes, the order of indication shall be $\theta_x\theta_y\theta_z$ based on the axial direction indicated in the dimensional drawing.)
- ⑦ Linearity
Linearity in closed loop control.
- ⑧ Capacitance
The capacitance of piezo-electric element is indicated as a representative value.
- ⑨ Stiffness
Rigidity in the travelling direction of the trolley table is indicated as a representative value.
- ⑩ Weight
Weight of the main unit of the Piezo Stage. (Cables and connectors are not included.)
- ⑪ Body material, surface treatment
Material and surface treatment method used for the Piezo Stage main unit.
Weight to carry in the accuracy inspection and measurement of the stage is prescribed in the internal standard, which is approximately 50 to 100 g as a reference. (This may slightly vary depending on items to measure.)
Please specify the actual load amount of your environment so we can make a measurement with the same load.

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
Calculation formula of the resonance frequency diagram

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Assume, for the spring mass model, that the rigidity of a piezo stage with an elastic hinge mechanism is K [N/m] and the mass of the moving section is M [kg], and that the resonance frequency f [Hz] when a load m [kg] is applied to this can be expressed in a formula of simple harmonic motion.

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Calculating the resonance frequencies f_1 and f_2 when loads m_1 and m_2 are applied in FEM analysis will find the apparent M and K values, and assigning these values in the above formula will lead to a simple harmonic motion formula for each piezo stage. Then, the resonance frequency f for an arbitrary value m will be found.

The figure of resonance frequency is a graphic presentation of this.

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The load applied in the FEM analysis is calculated with an assumed rectangular solid. Therefore, there may be a difference between the actual resonance frequency and the graph in the catalog depending on the shape or the center of gravity of the load or the position in which the product is installed on the table.

If the customer provides us with the information of the load, we will design an equivalent model and calculate the resonance frequency through FEM analysis.

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

Increase the voltage from 0 V to the rated voltage at a constant rate and then decrease the voltage from the rated voltage to 0 V at a constant rate. Repeat the above cycle 5 times and take the difference between the maximum displacement and minimum displacement in the hysteresis curve (piezo-electric hysteresis) of the 5th cycle as a stroke value for open loop.

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Input the minimum step position command for at least 5 steps successively in the positive (or negative) direction relative to an arbitrary point (a reference measurement point) and measure the actual position. The maximum travel distance divided by n (number of steps) is the resolution. The resolution shall be subject to inspection only for closed loop control with built-in displacement sensor.

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
* Resolution for open loop control (theoretical resolution)

Even loop control operations may achieve fine steps in nm by applying in small steps. However, this resolution represents the amount of fine movements caused by small applied voltage without reading of displacement sensor signals in the stage, so it is also referred to as "theoretical resolution." We define resolution as fine steps in a condition that the piezo stage is controlled by the controller by means of feedback control and it can actively control the position.

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
Repeat the positioning to an arbitrary point and measure the stop position at least 7 times and then obtain 1/2 of the maximum difference of the measurement values.

This measurement should be made at the center of the stroke and, as needed, at both ends of the stroke.

The repetitive positioning accuracy shall be subject to inspection only for closed loop control.

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
While checking the portion where the phase is changing by using an impedance analyzer, measure the frequency of the resonance point (a portion where impedance is the lowest).

There are also alternative measurement methods:

- Inputting rectangular wave and measuring the resonance frequency
- Inputting sine wave while frequency sweeping and measuring the resonance frequency (Measurement by impedance analyzer)

Piezo stage ▶ P.11
Piezo motor stage ▶ P.127
Force sensor ▶ P.153
Technical information
Characteristics of piezo-electric element
Mechanism and material of piezo stage
Capacitance type displacement sensor
Open loop and closed loop
Cable extension
Vacuum specifications
Installation direction of piezo stage
How to read the specification table of piezo stage
Calculation formula of the resonance frequency diagram
Operation mechanism features
Inspection method of piezo stage
Operating data example (Piezo stage)
PID gain adjustment example (Piezo stage)
Sample software (piezo stage)

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



Pitching : Perform positioning from the reference position (stroke end) to the maximum stroke in the same direction successively at certain intervals and use an autocollimator to measure the inclination in the vertical direction in each positioning point to obtain the maximum angular difference.

Yawing : Measure the inclination in the horizontal direction with the maximum stroke to obtain the maximum angular difference.

Rolling : Set an autocollimator at right angle to the stage travel direction to make a measurement.

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Perform positioning from the reference position (zero position) successively in one way at certain intervals and measure the position in the travelling direction.

Approximate the ideal line (linear function) by the method of least squares based on the measurement result, and take the ideal line as a reference to identify a point where the difference from the actual measurement becomes the highest.

The linearity shall be subject to inspection only for closed loop control with built-in displacement sensor.

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This data is an example.
 Positioning accuracy and dynamic characteristic performance depends on various conditions such as stage stroke, resonant frequency, weight and shape of the load, controller type, and PID gain adjustment.

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
Resolution · Positioning repeatability

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
Creep (drift) of piezoelectric elements in open loop operation

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


Sine wave motion in closed-loop control

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
Triangle wave motion in closed-loop control

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
This is an example of the change in stage response due to filter and PID gain adjustment.
 The effect of adjustments depends on various conditions such as the characteristics of the stage, weight and shape of the load, controller type, and command waveform.

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